





**UNIVERSAL THERAPY WRAP**  
INSTRUCTIONS FOR USE

## EASY TO USE.

- 1 Lock Hub securely into Wrap.



- 2 To turn on, simply hold button on the center of Hub for 2 seconds until a blue light is illuminated.




- 3 Press button once to begin low frequency vibration. Press button again to toggle between desired vibration levels.

Low frequency    High frequency    Pulse



- 4 Place Intellinetix Hub + Wrap on area in need, ensuring nodes of Hub surround targeted area of pain or muscle tension. Secure in place using hook and loop closures on Universal Wrap.

- 5 Apply vibration therapy for 25 – 45 minutes for maximum level of pain relief.



- 6 When finished, simply hold button down for two seconds to turn off.

- 7 Replenish battery by plugging USB charger into the charge port.

- 8 Use as needed to increase circulation, improve flexibility, aid recovery and maximize pain relief.

## EASY TO MAINTAIN.

Recharges in minutes with USB charger. Simply plug one end into AC wall adapter or powered USB port on alternative device and other end into the Hub. Lights will illuminate as unit charges. Blue light indicates full charge. Fully charged, Hub will operate up to 45 minutes on high frequency.

## EASY TO SERVICE.

Brownmed warrants this product to be free from manufacturing defects. If your Intellinetix Universal Therapy Wrap stops working within one year from date of purchase, we will replace them.

## EASY TO CARE FOR.

- Spot clean as needed.
- Lay flat and allow to air dry.
- Keep charger opening dry while cleaning.
- For adult use only.
- Do not dry clean.

### Need a replacement or extra parts?

Customer service is happy to help you.  
Please call us at **1-877-853-5518**.

Show us how **you** vibe.  
**#Intellinetix**

Learn more at [intellinetix.com](http://intellinetix.com)



**WARNING:** Lithium-ion Batteries and/or products that contain Lithium-ion Batteries can expose you to chemicals including cobalt lithium nickel oxide and nickel, which are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

**RM07548**

2018-0081 Rev. 1