

Toe Hold[®]



Adjustable Stabilization of First Metatarsal-Phalangeal Joint

- maintains the hallux in any corrective plane the physician desires (neutral, adducted, varus, dorsiflexed, or plantarflexed position) with appropriate placement of hook and loop strap
- improves hallux purchase during ambulation by stretching soft tissue contractures
- provides support for turf toe or broken hallux
- non-latex foam construction
- one piece design applies easily and stretches to fit without compromising stability

Effective Treatment Following

- first MPJ realignment
- hallus valgus procedure

Item #	Size
50040	S/M
50041	S/M Left
50042	M/L Right
50043	M/L Left

S/M for shoe widths AA – B
M/L for shoe widths C - EEEE

Material Content: Neoprene